


MASIMO

Radical-7™

Providing upgradable rainbow technology in multiple configurations for maximum clinical flexibility



Masimo Radical-7



- > Featuring “gold standard” Masimo SET® pulse oximetry, proven in more than 100 independent and objective studies to provide the most accurate and reliable SpO₂ readings during motion and low perfusion.
- > Upgradable Masimo Rainbow SET technology platform lets you add total hemoglobin (SpHb™) and total arterial oxygen content (SpOC™) as a factory-ordered option or through simple field-installed software upgrades.
- > Additional upgrades allow you to continuously and noninvasively measure carboxyhemoglobin (SpCO®), methemoglobin (SpMet®), and PVI™.



MASIMO RADICAL-7 AS A HANDHELD OR BEDSIDE CO-OXIMETER

The Radical-7 easily detaches from the bedside unit as a full-featured handheld pulse co-oximeter. The onboard rechargeable 4-hour battery and >10 day trending facilitate both hospital transport and spot check applications.



MASIMO SATSHARE


Upgrade your current system to Rainbow technology simply by linking the Radical-7 to a host monitor's oximeter input connector. Get Masimo SET SpO₂ readings on more than 100 multiparameter monitors and SpHb, PVI, SpCO and SpMet readings on the Radical-7 itself.



GRAVITY ACTIVATED ROTATION

Automatic display rotation (Gravity Activated) provides for vertical and horizontal positioning to maximize space utilization and visibility.

MASIMO RAINBOW SET MEASUREMENTS GIVE YOU A MORE COMPLETE PICTURE OF YOUR PATIENTS' PHYSIOLOGICAL STATUS




EASILY CONFIGURED TO MEET YOUR SPECIFIC CLINICAL NEED




NORMAL VIEW

Allows for continuous monitoring of all parameters along with a plethysmograph waveform and Signal IQ indicator (below pleth) to provide confidence in the quality of the values displayed



TREND VIEW

The quick-trend feature allows for one-touch access to vital parameter trending information to instantly evaluate patient condition and illness severity. Auto-scaling feature allows the y-axis to automatically change with fluctuations in the measurement value



SIMPLIFIED ALARM MANAGEMENT

One-touch menu provides quick access to the most commonly used features and allows for quick configuration and management of parameter alarm settings

- > The Radical-7 features a bright multi-color or monochromatic screen that is easy to read in either vertical or horizontal orientation—perfect for at-a-glance readings in a variety of clinical settings.
- > Patented 3D alarms provide thresholds for desaturations and perfusion index based on clinician-specified severity and frequency, enhancing patient safety.
- > Optional external wireless radio allows The Radical-7 to communicate with Masimo Patient SafetyNet, the remote monitoring and clinician-notification system that helps you keep at-risk patients safe on general care floors.



PERFORMANCE

MEASUREMENT RANGE

SpO ₂	0 – 100%
SpMet	0 – 99.9%
SpCO	0 – 99%
SpHb	0 – 25 g/dL
SpOC	0 – 35ml of O ₂ /dL of blood
Pulse Rate	25 – 240 bpm
Perfusion Index	0.02 – 20%
PVI	0 – 100%

OXYGEN SATURATION ACCURACY (%SpO₂)

Saturation	60 – 80%
No Motion	
Adults/Infants/Pediatrics	± 3%
Saturation	70 – 100%
No Motion	
Adults/Infants/Pediatrics	± 2%
Neonates	± 3%
Motion	
Adults/Infants/Pediatrics/Neonates	± 3%
Low Perfusion	
Adults/Infants/Pediatrics/Neonates	± 2%

PULSE RATE ACCURACY

Pulse Rate	25 – 240 bpm
No Motion	
Adults/Infants/Pediatrics/Neonates	± 3 bpm
Motion	
Adults/Infants/Pediatrics/Neonates	± 5 bpm
Low Perfusion	
Adults/Infants/Pediatrics/Neonates	± 3 bpm

CARBOXYHEMOGLOBIN SATURATION ACCURACY (%SpCO)*

Adults/Infants/Pediatrics	1 – 40% ± 3%
---------------------------------	--------------

METHEMOGLOBIN SATURATION ACCURACY (%SpMet)*

Adults/Infants/Pediatrics/Neonates	1 – 15% ± 1%
--	--------------

TOTAL HEMOGLOBIN ACCURACY (SpHb g/dL)

Adults/Infants/Pediatrics	8-17 g/dL ± 1 g/dL
---------------------------------	--------------------

RESOLUTION

Oxyhemoglobin Saturation (%SpO ₂)	1%
Carboxyhemoglobin Saturation (%SpCO), digital display	1%
Methemoglobin Saturation (%SpMet), digital display	0.1%
Total Hemoglobin (SpHb g/dL)	0.1 g/dL
Pulse Rate (bpm):	1 bpm

ELECTRICAL

STANDALONE

AC power requirements	100 – 240 VAC, 47 – 63 Hz
Power consumption	55 VA
Fuses	1 Amp, Fast Acting, Metric, (5x20mm), 250V

BATTERIES

HANDHELD

Type	NiMH
Capacity (battery life)	4 hours
Charging time	3 hours

STANDALONE (MODEL RDS-1B ONLY)

Type	NiMH
Capacity (battery life)	10 hours
Charging time	6 hours

ENVIRONMENTAL

Operating temperature	41°F to 104°F (5°C to 40°C)
Storage temperature	-40°F to 158°F (-40°C to + 70°C)
Operating humidity	5% to 95%, noncondensing
Operating altitude	500 mbar to 1060 mbar pressure -1000 ft to 18,000 ft (-304 m to 5,486 m)

PHYSICAL CHARACTERISTICS

Dimensions	8.9" x 3.5" x 2.1" (22.6 cm x 8.9 cm x 5.3 cm)
Standalone	3.5" x 10.5" x 7.7" (8.9 cm x 26.7 cm x 19.6 cm)

WEIGHT

Handheld	1.2 lbs (0.54 kg)
Docking Station .. (models RDS-1, 2, and 3)	2.5 lbs (1.14 kg)
(model RDS-1B)	3.1 lbs (1.41 kg)
Standalone	(models RDS-1, 2, and 3) .. 3.7 lbs (1.68 kg)
(model RDS-1B)	4.3 lbs (1.95 kg)

TRENDING

Provides 72 hours of trending at 2-second resolution for >10 days at 10 second resolution of SpO₂, SpHb, SpOC, SpMet, SpCO, Pulse Rate, and Perfusion Index with output to serial printer or other serial devices.

SpO₂ MODES

Averaging mode	2, 4, 8, 10, 12, 14, or 16 seconds
Sensitivity	Normal, APOD, and Maximum

ALARMS

Audible and visual alarms for high low saturation and pulse rate (SpO₂ range 1-99%, pulse rate range 30-235 bpm, SpCO 1-98%, SpMet 1-99.5%, SpHb 1-24.5 g/dL, PI .03-19%, PVI 1-99%)

DISPLAY/INDICATORS


Data display: %SpO ₂ , %SpMet, %SpCO, SpHb, pulse rate, pleth waveform, alarm status, trends, status messages, Signal IQ, perfusion index, MAX and APOD sensitivities, PVI and FastSat	
Display Color	Multi-color or monochrome blue
Type	TFT Color LCD
Pixels	480 x 272 dots
Dot Pitch	0.20 mm

OUTPUT INTERFACE

SatShare (RDS-1, RDS-1B); Serial RS-232 (RDS-1, RDS-1B, RDS-3); Nurse Call/Analog Output (RDS-1, RDS-1B, RDS-3); Philips Vuelink, Spacelabs Universal Flexport, (RDS-1, RDS-1B, RDS-3)


CLINICAL CONFIGURATIONS:

With a choice of docking stations for your Radical-7, you can select the connectivity configurations that work best for your clinical needs



RDS-1

Serial, analog, nurse call and SatShare connectivity. Optional extended battery provides battery life up to 10 hours. (RDS-1B)



RDS-2

Power Only.



RDS-3

Serial, analog and nurse call connectivity.

* SpO₂, SpCO, and SpMet accuracy was determined by testing healthy adult volunteers in the range of 60% - 100% SpO₂, 0% - 40% SpCO, and 0% - 15% SpMet against a laboratory CO-Oximeter. SpO₂ and SpMet accuracy was determined on 16 neonatal NICU patients ranging in age from 7 to 135 days old and weighing between 0.5 and 4.25 kgs. Seventy-nine (79) data samples were collected over a range of 70 - 100% SaO₂ and 0.5 - 2.5% HbMet with a resultant accuracy of 2.9% SpO₂ and 0.9% SpMet. Contact Masimo for testing specifications.