



Incorporated filter for viral and bacterial protection

Selector for adult or paediatric use

Main body in unbreakable techno polymer

Disposable jar

Ergonomic grip for ambidextrous use

Lightweight

B-Vac

Manual aspirator for secretions

B-Vac is an aspirator for secretions which is easy to use, portable, stable and compact and which generates a vacuum of -550 mmHg.

Just one hand is sufficient to generate the vacuum. An essential design allows for ambidextrous use.

In order to guarantee efficiency and precision the B-Vac incorporates a novel and practical system which permits selection for use on either adult or paediatric patients. In order to guarantee viral and bacterial protection the product has an easily replaceable filter. The 110 ° rotation of the tube allows for better access even in restricted areas. The use of the latest generation of techno polymers for the construction ensures a long life and easy sanitation. Further information on the comparison chart pg. 424.

B-Vac

RV71150B

S	CE	Comp 424	NEW	
E PS	P	10 x 40 mm	col 02	min 01
I interest	HC	LC	PD	PR