


Blue Splint Kit

Splints with flexible internal structure and Velcro® closure



A system created by Spencer to immobilize joints effectively, simply and practically. The Blue Splint immobilization system, made in five shapes, has changed its way of treating fractured or sprained joints. All splints can be modelled according to the necessity of each case and, once obtained the desired profile, the splint will offer extraordinary rigidity. To place the splint in a straight position it is sufficient to straighten by hand the inside structure. The Blue Splints are extremely light and could be overlapped with each other to create an infinity of configurations and solutions. They can be used for adults and children of every size. The versatility and the compactness that characterize the splints make them useful for the rescuer even in the most strange of situations. The arm and forearm splints are also used as an armrest during intravenous infusions. Perfect for ambulances, medical cars, mountain rescue, sport medicine, water rescue, orthopaedic centres, first aid and all kinds of emergency departments. The lightness of the complete kit and the reduced dimensions make it the most cumbersome and compact system on the market.

Technical details

Dimensions (bag) 640x100xh310 mm
Weight 1,5 kg

Latex free

JM80002A

BLUE SPLINT Kit

catalogue pg. 174

