

**SUPER BLUE**

Compact universal head immobilizer



Indispensable for head immobilization. The Super Blue can be used with spine boards and the scoop stretcher. It is made of high density foam, covered with water repellent material.

- Easy to use
- Washable and easy to stock
- The head supports have an anatomic profile and offer a stable and comfortable support
- The head support's position can be adjusted thanks to the use of a highly resistant type of Velcro®
- The head immobilizing straps are made out of polypropylene and are padded in correspondance to the part in contact with the patient's face
- The positioning of the base is very easy, thanks to the Velcro® fixed to the lower part of the base.
- The base is reinforced in the areas most exposed to use.

**TECHNICAL DATA**

Total weight (kg):	1 kg
Compatibility with spine board:	All
Fixation straps:	2
Type of fixation:	Velcro® / straps
Dimensions (mm):	384 x 21 x h290
Type:	Stuffed
Material:	Expanded PE / Spentex®
Manufacturing system:	Injection
Colour:	Blue
<b>Head blocks</b>	
Dimensions (mm):	245 x 80 x h160
Weight (kg):	0,5
Structure:	Stuffed pillows
Material:	Plastic compound
Manufacturing system:	Injection

Super Blue blue SH00110C

