

# EMERGENCY

Emergency bag empty with English writing



We have concentrated on the functionality of the product: a lightweight bag, easy to handle and easy to clean. Eliminating the superfluous, all you have left is the handiness of the product and the convenience of its price. Made of nylon, it has a central compartment where you will find on the inside a premodelled and semi-rigid compartment with three dividers to store a small oxygen tank and other resuscitation devices.

The closure of the top flap is fast and safe thanks to the Velcro® strap all around the front border.

There are also two side pockets with zipper.

The Emergency stands out for the global character of its range: it is available in fact in different languages.

## Technical details

Dimensions	525X215Xh230 mm
Capacity	26 L
Weight	0,7 kg
Material	Ripstop 600 D
Carry strap	Yes
Lateral pockets	Yes
Rigid double bottom	Yes
Bottom reinforcement	Yes
Reflective insert	Yes
Int. material	Alveolar polypropylene
Personalization	Yes, hot print 180x80 mm

## Latex free

CB08001A

EMERGENCY BAG