

Green
spec-2™

The Most Extensive Range of Fibre Optic Laryngoscope Blades with FibreClip™ Replaceable Light Guide



- Widest choice of Macintosh and Miller Blades and the unique Oxyport paediatric blades
- Incomparable cold, bright light - 2800 LUX
- Simple, trouble-free use. No screws, no lost parts
- Autoclavable FibreClip™ light guide, tested to withstand over 500 cycles
- Stainless steel parts guaranteed for 5 years

Conforms to ISO 7376:2009
Approvals: ISO 9001:2008, ISO 13485:2003
FDA registered manufacturer



Technical Specifications

Blades				
Size	Use	Overall Length (mm)	Distal Width (mm)	Order No.
Mac0	Neo-natal	77	9	4230
Mac1	Neo-natal	90	9	4231
Mac2	Infant	100	12	4232
Mac3	Adult	130	15	4233
Mac4	Adult	155	15	4234
Mac5	Adult	172	15	4235

Size	Use	Overall Length (mm)	Distal Width (mm)	Order No.
Miller00	Neo-natal	64	11.5	4249
Miller0	Neo-natal	75	11.5	4250
Miller1	Infant	102	11.5	4251
Miller2	Child	155	14	4252
Miller3	Adult	195	14	4253
Miller4	Adult	205	14	4254

Size	Use	Overall Length (mm)	Distal Width (mm)	Order No.
OxyMiller00	Neo-natal	76	11.5	4240
OxyMiller0	Neo-natal	87	11.5	4241

Sets Include*		
Green SpeX™ Handle	GreenSpec 2™ Blades with FibreClip™ bundle	Order No.
Medium size	Mac 1,2 & 3	4260/X
Medium size	Mac 1,2,3 & 4	4261/X
Medium size	Mac 2,3 & 4	4262/X
Medium size	Miller 00,1, 2 & 3	4263/X

* Additional combinations available upon request

FibreClip™ Light Guide			
Size	Use	Overall Length (mm)	Order No. Pack of 2 FibreClip™ Light Guides
Mac0	Neo-natal	55	423002
Mac1	Neo-natal	55	423102
Mac2	Infant	67	423202
Mac3	Adult	91	423302
Mac4	Adult	117	423402
Mac5	Adult	117	423502
Miller00	Neo-natal	37	424902
Miller0	Neo-natal	43	425002
Miller1	Infant	71	425102
Miller2	Child	106	425202
Miller3	Adult	147	425302
Miller4	Adult	157	425402

Distributed by

Product	Description	Ordering No.
Tru-LED™ Adult Handle set	Tru-LED™ 32 mm diameter ergonomical handle with rechargeable battery cartridge and power supply in clear plastic case	4220 E,A,B or U/J*
Tru-LED™ Infant Handle set	Tru-LED™ 28 mm diameter ergonomical handle with rechargeable battery cartridge and power supply in clear plastic case	4221 E,A,B or U/J*

Sets Include		
Tru-LED™ Handle	GreenSpec2™ Blades with FibreClip™ bundle	Order No.
Adult Size	Mac 1,2 & 3	4260/TL- E,A,B,U
Adult Size	Mac 1,2,3 & 4	4261/TL- E,A,B,U
Adult Size	Mac 2,3 & 4	4262/TL- E,A,B,U
Adult Size	Miller 00,1,2 & 3	4263/TL- E,A,B,U
Infant size	Miller Oxy 00,0 & 1	4265/TI - E,A,B,U

* Legend: E = EU plug, A = Australian plug
B = British plug, U/J = USA/Japanese plug

GreenSpeX™ Handles			
Size	Finish	Batteries	Code
Medium	Knurled	C	4211X
Mini	Knurled	AA	4212X
Stubby	Knurled	AA	4213X

Bulbs		Order No. Box of 6
Green SpeX™ Handle	Xenon 2.5V 0.7A 4.8 O.D. M3.5x0.35	512501FO/6

From Truphatek®'s Warranty Statement:

All Truphatek® products are unconditionally guaranteed to be of good quality and free from manufacturing defects in material and workmanship for a period of 5 years from date of manufacture. Replacement product will be provided without charge for any Truphatek® product not meeting this warranty.

FibreClip™ replaceable fibre-optic light guides are guaranteed to perform for a period of 18 months from date of manufacture"

Full text of Truphatek®'s warranty can be viewed at our website www.truphatek.com or by request from Customer Service Desk at helpdesk@truphatek.com

Green Spec2® is a patented product. Green SpeX™, Green Spec2® and Tru-LED™ are trademarks of Truphatek® and should not be used without prior written permission from Truphatek International Ltd.



Truphatek International Ltd.
POB 8051 Netanya 42504, Israel
Tel: 972-9-8851155, Fax: 972-9-8851212,
helpdesk@truphatek.com