


SPENCER 122

Four fold stretcher



Thanks to reduced clutter, the stretcher is used as an emergency stretcher and also in conditions where storage space lacks.

- Complete with two orange quick-release restraints
- Reflecting mark to increase nocturnal visibility
- It is also available in the decontamination version

Technical details

Foldable stretcher	Yes
Lenght open	2010 mm
Lenght folded	1000 mm
Width open	400 mm
Width folded	105 mm
Height open	140 mm
Height folded	165 mm
Weight	6.5 kg
Load capacity	170 kg
N. of carrying points	4
Feet	4
Handles	Yes
Frame	Aluminium
Reflective mark	Yes
Steel folding mechan.	Stainless
Two belts	Yes
Sheet color	Orange
Sheet material	Spentex®
Electronically welded	Yes
Suitabòe for decontam.	ST00126A

Latex free

ST00122A

Spencer 122 orange

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
E-mail: info@spencer.it Web: <http://www.spencer.it>

This document has been automatically generated by www.spencer.it on 2010/01/14 12:12:19, and is referred to the last production lot. The informations furnished have been last updated on 2009/12/01. To verify the correctness of this document, please contact our Quality Dept. at quality@spencer.it