

# Porta Vac

Compact suction device



The impressive simplicity of the shapes of this model, together with the distinguished technical characteristics which guarantee extreme functionality, immediately attract attention.

The rectangular vacuumeter perfectly integrates with the shapes of the body and is placed in a protected position, perfectly visible from above.

The high-capacity engine does not need any maintenance and easily creates high grades of vacuum (-550 mm Hg) and high capacity (23 L/min).

The Porta Vac is equipped with a local power source plug, allowing quick recharging of the batteries and/or functioning with net-energy.

Porta Vac can be supplied with a comfortable carrying bag, which can contain not only the suction device, but all accessories as well.

Another important accessory is the metal wall-support for fixation on emergency vehicles.

## Technical details

Dimensions	370x100x240 mm
Weight	5,5 kg
Max Vacuum	- 550 mm Hg
Max capacity	23 L/min
Power	12 / 110 / 220 V
Container	Yes
Filter	1000
Silicone tube	Yes
Use temp	0°C to +40°C
Storage temp	-10°C to +50°C
Humidity	95 %
Wall support	Yes
10G support	Yes
Light signals	Yes

Latex free

SC60001A

PORTA VAC

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
E-mail: [info@spencer.it](mailto:info@spencer.it) Web: <http://www.spencer.it>

This document has been automatically generated by [www.spencer.it](http://www.spencer.it) on 2010/04/23 12:12:19, and is referred to the last production lot. The informations furnished have been last updated on 2009/12/01. To verify the correctness of this document, please contact our Quality Dept. at [quality@spencer.it](mailto:quality@spencer.it)