

MAXIZAINO

Professional Emergency Rucksack



The inside of the rucksack is divided into two rigid compartments and the back is covered with two large bands of velcro for the pouch to attach onto. The inside of the front flap has a big pocket divided into 4 zip closed compartments. The bottom area is covered with an expanding pouch to optimise the space. On the outside of the front flap there are two big pockets with zip closure. The handle on the top and on the side make the rucksack easy to carry as a bag.

TECHNICAL DATA

Dimensions	500 x 220 x h450 mm
Volume	50 L
Weight empty	2.3 kg
Weight with contents	11.9 kg
External fabric	Nylon Oxford
Back fabric	Nylon Oxford
Internal fabric	Polyester
Structural material	Corrugated PE
Inner fabric	Polyester
Structural material	Corrugated PE
Handle cover	TPU
Strengthened base	Yes
Padded shoulders	Yes
Frontal compartments	Yes
Inner pockets	2
Gas cylinder fixing	Yes
Zip pulls	Yes
Reflective inserts	Yes
Maxizaino	CB04040