

CROSS

Roll-in stretcher trolley with adjustable heights black board



Made in lightweight stainless steel and high density polyethylene, the Cross is practical as well as comfortable to use. Bigger and stronger wheels make the difference in moving patients on rough terrain, as stability and balance are enhanced. The side bars are safely locked thanks to a quick-release mechanism, which prevents accidental unlocking. One of the most important features of the Cross is its flexibility: the agile stretcher can be adapted to any type of vehicle.

Lightweight, fast, precise, multi-skilled and consistent to name a few are the wonderful characteristics of this revolutionary stretcher, which is the result of the experience of Spencer's research and development efforts.

It is the attention to detail that makes a stretcher special. Discover how ergonomically detailed the Cross is.

IV Pole holder	yes
Wheels Ø 200 mm	2 front fixed, 2 rear auto directional turning 360° with break
Patient loading board	moulded rigid monoshell in high density polyethylene
Command labels	yes
Cardio-pulmonary resuscitation	yes
Adjustable head/backrest	yes
Opening angle backrest	min. 75° (and at least 5 possible positions)
Trendelenburg board	adjustable height at least 900 mm
Trendelenburg board	possibility to raise at least 15°

Technical details

Length	1970 mm
Shortened length	1465 mm
Width	570 mm
Loading height	from 480 to 750 mm
Weight	38,7 kg
Load capacity	160 kg
Folding side bars	yes
Length	630 mm
Height	180 mm

CR10054E

CROSS

