

URRA' SPLINT

Semi-rigid splints size L



This particular immobilization system has exceptional characteristics.

A thin layer of aluminium is covered with expanded closed cell polypropylene, impermeable to fluids and blood, and it is X-ray compatible.

When it is longitudinally nerved, the splint becomes a semi-rigid support for the distal parts of the limbs.

The surface is soft and sufficiently resistant to be sanitised and re-used over again.

Thanks to its small size, it can be placed inside all types of first aid kits.

Available in different sizes.

Technical details

Lenght	905 mm
Width	5 mm
Height	108 mm
Material	Aluminium / polypropylene
Use temperature	From -10°C to +55°C
Storage temperature	From -20°C to +60°C

Latex free

JM90000A

URRA' SPLINT
CM.90X11