

Spencer SX

Scoop stretcher 3 blades with 3 belts



SPENCER SX is a device which enables the lifting and loading for a patient onto immobilising or transport devices such as spinal boards, vacuum mattresses or transport stretchers.

It is essential for the loading trauma patients by two rescuers: with use of this device patient can be picked up without changing position and the trunk and limbs can maintain the supine position while loading.

The structure and mechanism of the SPENCER SX offer maximum reduction of the movement required. The SPENCER SX does not offer any type of isolation. No action on behalf of the patient is expected.

This model is supplied with three belts.

Technical details

Material	Aluminium
Weight	9.9 Kg
Load capacity	170 kg
Height	65 mm
Width	435 mm
Minimum length	1665 mm
Maximum length	2168 mm
Folded length	1200 mm
Finishing	Anodized

Latex free