

# B-bak Pin

## Spine board complete of pins



A special B-bak, very much appreciated by rescue teams, it puts together its famous technical characteristics and a built in system that allows to fix the restraints to the board. The side handles are translucent incorporated elements, they allow you easy buckling-unbuckling operations for a fast and safe anchorage of the patient's restraints with shackle.

High density completely recyclable polyethylene, moulded into a unique, one piece shell; this technique guarantees the - B-bak the following exclusive characteristics:

- rigidity and lightness
- resistance to bumps and corrosion
- completely waterproof (does not absorb liquids)
- immune to infiltrations
- easy to clean (water and soap are sufficient)
- totally X-ray compatible (the patient can be removed from - - the spine board even after carrying out X-rays or CAT scan investigations).

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
E-mail: info@spencer.it Web: <http://www.spencer.it>

This document has been automatically generated by [www.spencer.it](http://www.spencer.it) on 2010/01/14 12:12:19, and is referred to the last production lot. The informations furnished have been last updated on 2009/12/01. To verify the correctness of this document, please contact our Quality Dept. at [quality@spencer.it](mailto:quality@spencer.it)

### Technical details

Dimensions	410x50xh1840 mm
Material	HDPE
Weight	8 kg
Load capacity	180 kg
Temp. of use	-20 +60 C°
No. of handles	14
Flexybility with 180kg at -10 C°	25 mm
Flexybility with 180kg at 20 C°	40 mm
Flexybility with 180kg at 50 C°	42 mm

X-ray compatible  
Watertight  
Latex free

ST02061B

catalogue pg. 122



B-bak pin