

SKID-OK

Lightweight evacuation chair



The most compact and technologically advanced evacuation chair by Spencer

Skid-Ok reaches high quality standards through the use of innovative materials and components that make it lightweight and manoeuvrable, able to operate without compromise in any operational scenario.

It features an aluminium frame, welded and polished and an adjustable seat, reinforced with Spentex®.

The innovative self braking system on the slides allows the operator a perfect control during the descent of the stairs.

The Skid-Ok chair is supplied with two polypropylene 50 mm straps.

TECHNICAL DATA

Width	520 mm
Length	900 mm
Height with fully extended lift bar	1540 mm
Height closed	1040 mm
Width closed	175 mm
Weight	10 kg
Loading weight	150 kg
Surface treatment	Polished
Front wheels	n°2 Ø 100 x 32 mm
Rear wheels	n°2 Ø 150 x 30 mm
Sheet	Spentex®
Sheet colour	Black/Red
Frame	Aluminium
Frame colour	Chromed
Belts	50mm Polypropylene

Latex free

Skid-Ok SK20003E

