


HIRVIN

Transfer sheet, high resistant double layered nylon with 6 handles



For excellent, quality, hard work and love for details are necessary. With this craftsmanship we are able to achieve the perfection that distinguishes our Hirvin injury transportation sheets.

Scrupulous controls guarantee that every sheet responds to high standards of quality. In fact the Spencer logo is a worldwide symbol of high quality.

During the realisation of the Hirvin injury transportation sheets a lot of attention was given mainly to the ergonomic aspect, considering the actual structural dynamic conditions as during use by the operator.

The invisible details that express it's robustness and manageability.

The composition, weight and framework of the material:

- 550 D nylon coated twice with poly-vinyl (Spentex®)

The composition, weight and finishing of the tape:

- 50 mm nylon tape with breaking load in traction 600 kg
- finishing to improve softness Types and configuration of the reinforcements:

- cross lock stitching that exploits the new electric control system of the tension of the thread (active tension), the length of the point and the width of the crosses.

The safety stitching:

- doubled in correspondence to all joining points subject to traction; all carried out with nylon fibre thread.

The quality of the thread:

- the thread, as with the stitching system it is very important, because the stitching is made to last; a resistant, unbreakable nylon thread is used

Technical details

Dimensions	1850 X 760 mm
Weight	1 Kg
Material	Nylon & Polyvinyl
N° handles	6 (3 each side)
Load capacity	150 kg

QM40100A

HIRVIN