



MR conditional laryngoscope system



Intubate safely within the MR environment

- ▲ **MRI Conditional** – Tested by Independent Laboratory for up to 3 tesla
FDA approved: D138670
- ▲ **Gold plated FiberClip™ and black handle** for easy identification and to avoid fatal mistakes.
- ▲ **Competitive price** to avoid hospital buying resistance.
- ▲ **Patented system**
 - Conventional blade with fiber optic guide.
 - Unique inexpensive battery packs.
- ▲ **Modular sets** - Numerous options of sets customized according to clients' needs.

Tested by:

SHELLOCK
R & D SERVICES, INC.



Tru-MR™ products have been tested to be MR- Conditional by Shellock R&D Services, Inc., (www.magneticresonancesafetytesting.com) for magnetic field interactions at 3-Tesla with highest spatial gradient field of 720- Gauss/cm or less.

MRI Conditional up to 3 Tesla
Conforms to ISO 7376:2009
Approvals: ISO 9001:2008, ISO 13485:2003
FDA registered manufacturer



TRUPHATEK[®]
www.truphatekchina.com

*“Complex airway management (e.g. fiber optic intubation) should be performed in a controlled environment outside of zone IV. **Alternative MRI-safe/conditional airway devices should be immediately available in the MRI Suite**”*

Practice Advisory on Anesthetic Care for Magnetic Resonance imaging, Anesthesiology 2009; 110: 459-79

Tru-MR™ Technical Specs:

Tru-MR™ Sets	
Code	Description of content
5442M	One handle, 5 blades (Miller 0 and 1, Mac 2, 3 and 4) and 4 Battery Packs
5443M	One handle, 6 blades (Miller 1, 2 and 3, Mac 2, 3 and 4) and 4 Battery Packs
	One handle, 1 blade and 2 battery packs (All blade sizes Macs and Millers available)

Tru-MR™ FiberClips™	
Code	Product Description
524102M - 524402M	Mac 1- 4, Pack of 2
525002M - 525402M	Miller 0 - 4, Pack of 2
Tru-MR™ Accessories	
Code	Product Description
53102M	2 Battery Packs



Distributed by:



Truphatek's **warranty statement** can be viewed on our website: www.truphatek.com
For further information please contact: info@truphatek.com

Tru-MR, FiberClips and Truphatek are trademarks of **Truphatek International Ltd.** and should not be used without prior written permission.

MR-2-05/14 Rev1

