

NRT

Neil Robertson type



A rescue stretcher specifically designed for the recovery of victims in difficult and cramped situations, especially where the patient has to be raised and carried to safety. Made with semi-rigid splints and high resistant cotton. It guarantees utmost safety for the patient also during lifting.

Thanks to a robust construction and a practical storage system it becomes indispensable for rescue in:

- mines
- industries
- building sites
- ships

TECHNICAL DATA

Dimensions: 1080 x 20 x h1830 mm

Weight: 7 kg

Load capacity: 220 kg

NRT white

ST30150A

