

Vacu Aide Compact

Portable compact suction unit



Vacu Aide Compact is the expression of a high quality technology, studied, projected and produced to face all emergency situations with efficiency.

You will be surprised with its compactness, lightness and versatility. The engine which does not require any maintenance, guarantees a capacity of 29 L/min a -600 mm Hg.

Thanks to the Quiet-Drive technology, if used for home care or hospitals you will appreciate its silence. All noise and vibrations are minimal.

Vacu Aide Compact is a new technological development: we have reduced the power consumption to guarantee more than one hour of use with battery supply and have improved performance to assure a whole day of continuous function with a 12 Vdc or a 220 V supply.

Its professional aspect will naturally fit into ambulances, backpacks or hospitals.

Great performance within one, extremely light and compact portable suction pump.

Great performance that other competitors will have to deal with.

Technical details

Dimensions (mm)	229x180x205
Weight	3.1 kg
Max Vacuum (mm Hg)	Adjustable from 80 a 550
Max capacity	27 ± 2 L/min
Source	Battery
	12 volt
	220 volt
Autoclavable container	Si
Vol. supplied container	800 ml
Filter	Si
PVC tube	Si
Temp. of use	0°C + 40°C
Temp. for storage	-4°C + 70°C
Humidity	95%
Wall support	Yes
10G Support	Yes
Light signals	Yes
Colours available	Grey / Yellow

SC73010B	Vacu Aide Compact grey
SC73100B	Vacu Aide Compact yellow