

## EMSOFT

### Vacuum splints



The EMSOFT splints are made of coupled vinyl with an ultrasonic welding process. They are easy and fast to model. They can also be used for exposed fractures. The standard version has Velcro® closure straps.

X-rays can be performed without removing the splints.

The valve, welded onto the PVC sheet, guarantees a strong hold; thanks to a one way system, it avoids air from entering when the vacuum pump is removed.

### TECHNICAL DATA

Dimensions EMSOFT Kit three sizes:

910 x 170 x h650 mm

Weight EMSOFT Kit two sizes: 2,8 kg

Weight EMSOFT Kit three sizes: 3,6 kg

Functioning temperature: from -5 to +40 °C

Storage temperature: from -10 to +50 °C

The EMSOFT kit is complete with:

- splint size S (750 x 600 mm)
- splint size M (900 x 700 mm),  
only into kit including three sizes
- splint size L (1100 x 800 mm)
- vacuum pump
- Ripstop case

Weight (kg) size S: 0,6

Weight (kg) size M: 0,9

Weight (kg) size L: 1

External covering: Vinyl

Internal material: Polystyrene

HFV: Yes

X-ray compatible: Yes

Latex free: Yes

Temperature of use (°C): -5 +40

Temperature of storage (°C): -10 +50

EMSOFT Kit three sizes

QM22707A

