



NO.1 CHOICE OF HOSPITALS

Vario

SUCTION PUMP



Flexible and mobile

- Lightweight and compact
- Ultra-quiet and low vibration
- Available in different versions
- Rechargeable-battery or electric operation

Main applications:
Emergency department
Hospital ward
Thorax drainage

Vario

SUCTION PUMP

The Vario may be used in many different situations – on the hospital ward, in the emergency department, in the postoperative phase, in the nursing home, in the doctor's practice or at home. It is suitable for many different needs.

■ Latest technology

The QuatroFlex™ is a unique piston/cylinder system made from plastic, which operates with four vacuum-creating modules. This creates constant suction, and the Vario produces very little noise and vibration. The membrane vacuum regulator with safety knob guarantees a rapid build-up of vacuum and allows accurate regulation of the vacuum.

















■ Flexible and mobile

For use on the hospital ward, the Vario can be attached to a mobile stand or onto standard rails. It can also be easily carried, thanks to its light weight and handle.

■ Versions

Medela offers the Vario as a high- or low-vacuum pump, with or without rechargeable battery.

Technical data

	Vario 8	Vario 8 AC/DC	Vario 8 c/i AC/DC	Vario 18	Vario 18 AC/DC	Vario 18 c/i AC/DC
 vacuum	low vacuum - 90 mbar - 9 kPa - 68 mmHg	low vacuum - 90 mbar - 9 kPa - 68 mmHg	low vacuum - 90 mbar - 9 kPa - 68 mmHg	high vacuum - 75 kPa - 563 mmHg	high vacuum - 75 kPa - 563 mmHg	medium vacuum - 55 kPa - 413 mmHg
 flow	low flow 8 l/min.	low flow 8 l/min.	low flow 8 l/min.	low flow 18 l/min.	low flow 18 l/min.	low flow 18 l/min.
AC versions						
Vario 8 and Vario 18	AC	230–240 V	50 Hz	90 W		T 0.8 A 5 x 20 mm
Vario 8 and Vario 18		230–240 V	60 Hz	90 W		T 0.8 A 6.3 x 32 mm
Vario 8 and Vario 18		120 V	60 Hz	70 W		T 1.25 A 6.3 x 32 mm
AC/DC versions						
Vario 8 and Vario 18	AC/DC	100–240 V	50/60 Hz	80 VA		T 1.0 A 5 x 20 mm
Vario 8 c/i and Vario 18 c/i	AC/DC					
 Operation	+5 ... +40 °C	 Transport/Storage	+5 ... +50 °C	P	70 ... 106 kPa	
 Dimensions	380 x 170 x 285 mm			RH	20 ... 95 %	
 Weight	AC 3.5 kg AC/DC 4.2 kg			CLASS	 AC	 AC/DC
					ISO 9001, ISO 13485	
					MDD (93/42/EEC) Class IIa 0123	
					CE	



Vario 18 (AC/DC) version

© Medela AG/177.0077/02.06/B

Contact us – we will be pleased to assist you.

Switzerland & International Sales

Medela AG, Medical Technology
Lättichstrasse 4b, 6341 Baar
Phone +41 (0)41 769 51 51, Fax +41 (0)41 769 51 00
e-mail: info@medela.ch, www.medela.ch